

Sakura

# Tissue-Tek® NEUTRA-GUARD™

Formaldehyde Safety Products for the Operating Room



**NEUTRALIZE THE PROBLEMS  
OF FORMALIN DISPOSAL AND CLEANUP**



**SAKURA**

Sakura

# Tissue-Tek® NEUTRA-GUARD™

## Formaldehyde Safety Products

### Three unique products help protect your health in the OR

- Control formalin spills in the operating suite
- Reduce personnel exposure to dangerous formalin vapors
- Certified technology for treating formalin waste
- Complies with federal, state, and local ordinances

#### Neutra-Pads™

- Ideal absorbent pads to place under formalin jars
- Neutralize formalin and reduce noxious odors
- Special pad backing prevents formalin leaks



Product Code	Quantity/Pad Size
4295	100 Pads 20.3 cm x 25.4 cm (8 in x 10 in)

#### Neutra-Wipes™

- Make it easy to wipe up and neutralize small formalin spills
- Convenient ready-to-use pop up packaging
- Perfect liner for tissue transporting trays, storage areas



Product Code	Quantity/Wipe Size
4278	Starter Pack 1 Dispenser, 1 Roll (60 wipes per roll) 30.5 cm x 31.8 cm (12 in x 12.5 in)
4279	4 Rolls

#### Neutra-Form™

- Effectively controls and manages large formalin spills
- The spill area is contained and neutralized
- Rapid gel action facilitates safe disposal



Product Code	Quantity
4304	2 Containers (each container neutralizes larger spills of 10% formalin or 2% glutaraldehyde)



*Proven Reliability*

Visit our web site at [www.sakuraus.com](http://www.sakuraus.com)

Sakura Finetek U.S.A., Inc.  
1750 West 214th Street  
Torrance, CA 90501 U.S.A.  
Phone: (800) 725-8723