

# Introducing another Superior Station



 **Tissue-Tek<sup>®</sup> Accu-Edge<sup>®</sup>**  
SAKURA  
Grossing Station

# Accu-Edge® PathPRO™ Grossing Station

Sakura offers a complete line of products for sample preparation and anatomical grossing. Accu-Edge® Grossing Stations provide the opportunity to customize a grossing station to fit a lab's specific needs. Sakura's grossing stations are built from surgical grade stainless steel and provide standardized performance while enhancing grossing capabilities.

The Accu-Edge® PathPRO™ is another superior grossing station. It is equipped with both downdraft and backdraft exhaust systems providing greater ventilation for safety.

The Accu-Edge PathPRO Grossing Station is available as an adjustable or fixed height unit. The elevating model #0137 has the ability to elevate from 32" to 44" to accommodate all lab personnel. Model #0138 is positioned at a standard counter height of 36".

## Standard Features

- Choice of fixed or adjustable height settings
- 2 ventilation systems
- High-quality stainless steel construction
- Vacuum breaker-protected water supply
- ½ horsepower disposal
- Polyethylene dissecting board
- Magnetic instrument holder
- Magnifier light, deck mount
- Paper towel dispenser
- Spray hose assembly
- Shelving



## Accu-Edge® Grossing Tools

Accu-Edge® grossing tools include specially designed grossing boards available in three sizes—standard, large, and wells—to accommodate a wide range of tissue types. Featuring two adjustable wells that ensure grossing to an exact thickness, these boards enable complete standardization of the grossing process.



In addition, Accu-Edge trimming knives and scalpel blades provide the sharp cuts necessary for exact grossing techniques. Tissue tampers allow a firm hold on specimens while making horizontal cuts.

The Accu-Edge product line includes three grossing forks designed to hold specimens including skin, arterial, or colon samples, and any specimens with layers, while maintaining perfect specimen form and integrity.



Trusted Innovation™

**Sakura Finetek USA., Inc.**  
1750 West 214th Street  
Torrance, CA 90501 U.S.A  
Phone: (800) 725-8723