

 **Tissue-Tek® Accu-Edge®**

SAKURA


Grossing Station



Trusted Innovation™

The most trusted name
in pathology

DESIGNED WITH YOU IN MIND



Instrument Holder




GFIC Waterproof Receptacles



Large sink



Base Leg Assembly



Dimensions: 83" (h) x 48" or 60" (w) x 29" (d)

OPTIONAL ACCESSORIES

- Air and Gas Valves
- Base Leg Assembly with Cabinet
- Camera, Video or Dictation Mount Facility
- Digital or Economy Scale
- Disposal, 1/2 HP, 3/4 HP or 1 HP
- Dissecting Area Rinse
- Eyewash Station
- Filter, Activated Carbon or Potassium Permanganate
- Microphone on Flex Arm
- Monitor and Keyboard Stand
- Self Contained Ventilation
- Video Camera Arm
- Voice Activated Dictation
- Variable Length either 48" or 60"

PRODUCT CODE:

- Pathology Workstation- 0115
- Base Leg Assembly with Cabinet- 0121



Sakura Finetek USA, Inc.
1750 West 214th Street
Torrance, CA 90501 U.S.A
Phone: (800) 725-7823
www.sakuraus.com