



# Tissue-Tek® Sectioning System

## Cutting-Edge Partners Block after Block Ribbon After Ribbon

Together, **Accu-Edge® Disposable Blades** and the **Accu-Cut® SRM™ 200 Rotary Microtome** create the ultimate in "system compatibility." Expect outstanding performance - section after section, ribbon after ribbon, time after time.

### Accu-Edge® Disposable Blades

- Manufactured from the highest quality materials
- Platinum and chromium coatings reinforce blade edge strength and durability
- Unique coating process enables smooth, precise, consistent sections
- Ultrasharp, uniformly consistent, undisputed in quality



### Accu-Cut® SRM™ 200 Rotary Microtome

- Stable, durable instrument
- Smooth-running handwheel
- Three-axis specimen orientation
- With or without specimen retraction
- Specimen trimming feature
- Lateral displacement of blade holder base
- Minimal maintenance required
- Accessories included for use with low and high-profile blades



### Accu-Cut® SRM™ 200 Rotary Microtome

#### Dimensions and Weight

Height:	29.5 cm (12 in)
Width:	40 cm (16 in)
Depth:	47 cm (18.5 in)
Weight:	29 kg (64 lbs)

**Thickness Setting Range** 0.5 - 60  $\mu$ m

### Accu-Cut® SRM™ 200 Microtome Code

(Standard accessories included with system)

with specimen retraction	1429
without specimen retraction	1429N
Blade Holder Base Assembly	1434
Holder for Disposable Blades	1435
High-Profile Back-Plate	1436
Low-Profile Back-Plate	1437
Orienting Clamp Attachment	1438
Universal Cassette Clamp	1439
Maintenance Kit	1440
(Microtome oil included)	

### Accu-Edge® Disposable Blades Code

Accu-Edge® High-Profile Blades	4685
Accu-Edge® Low-Profile Blades	4689
Tissue-Tek® Paraform®	
Low-Profile Blades	7089
Tissue-Tek® Feather®	
Low Profile Blade Dispenser	4980
Disposable blades/case included	



Trusted Innovation™

Sakura Finetek USA, Inc.  
1750 214<sup>th</sup> Street  
Torrance, CA 90501  
Phone: (800) 725-8723  
[www.sakuraus.com](http://www.sakuraus.com)